1 Introduction

Mtzdump is a CCP4 "jiffy" program for extracting MTZ header information into a human readable form. This class is designed to provide this functionality at the Python interface level.

2 Use Cases

2.1 Simple: Display the Contents

This will need to be able to:

- Display the contents of an MTZ file.
- Raise appropriate errors if the input file is not MTZ.
- Ascertain the "kind" of mtz file.

Mtzdump has output like:

```
* Title:
* Base dataset:
        HKL base
* Number of Datasets = 1
* Dataset ID, project/crystal/dataset names, cell dimensions, wavelength:
      1 Unspecified
        12287_1_E1
        Unspecified
51.6515 51.6515 157.6742 90.0000 90.0000 90.0000
* Number of Columns = 18
* Number of Reflections = 6199
* Missing value set to NaN in input mtz file
* Number of Batches = 10
* HISTORY for current MTZ file :
From MOSFLM run on 5/6/06
* Column Labels :
H K L M/ISYM BATCH I SIGI IPR SIGIPR FRACTIONCALC XDET YDET ROT WIDTH LP MPART FLAG BGPKRATIOS
* Column Types :
HHHYBJQJQRRRRRRIIR
* Associated datasets :
000000000000000000
* Cell Dimensions : (obsolete - refer to dataset cell dimensions above)
  51.6515 51.6515 157.6742 90.0000 90.0000 90.0000
* Resolution Range :
```

```
0.00036 0.11111 ( 52.559 - 3.000 A )

* Sort Order :

0 0 0 0 0 0

* Space group = 'P43212' (number 96)

Batch number:

1
Batch number:
2
.... etc ....
```

then a small amount of "overall" information from the whole of the reflection file. This is run as "mtzdump hklin foo.mtz".

3 API

3.1 Output Dictionary

The output dictionary will contain:

- column_labels the column labels.
- column_types the (single character) column types.
- spacegroup the spacegroup.
- datasets a list of pname/xname/dname tokens.
- dataset_info a dictionary of unit cell and wavelengths associated with different datasets.

So far the output dictionary looks like: